

hearX STK
self test kit

The ideal teleaudiology solution

Effortless and precise remote hearing support




The hearX Self Test Kit (STK) is a tablet-based hearing solution that enables contactless hearing testing.

Make informed treatment suggestions and support patients with interventions like hearing aids, with this low- or no-touch hearing assessment solution, ensuring patient safety without compromising diagnostic accuracy.

Customize the kit with additional features tailored to your specific needs, enhancing the flexibility and functionality of this hearing assessment solution.

Expand your healthcare reach with these use cases



Home testing

No-touch, quick and user-friendly at-home hearing assessment.



Curb-side testing

Offers a low-touch, drive-through hearing assessment service.



Counter-side testing

Low-touch, quick and easy in-store hearing assessment.

Customize your Self Test Kit

Ideal for low- or no-touch alternative audiology care that includes:



Pure tone audiometry

hearTest™ is a CE-marked, FDA registered pure tone audiometer, used with calibrated headphones, for rapid automated audiometry.



Speech-in-noise

hearSpeech measures difficulty understanding speech-in-noise using a binaural digits-in-noise test. MCL/UCL testing was included to determine the dynamic reach of patients for hearing aid fitment. Access the [clinical study](#) to find out more.



Ear risk assessment

hearRisk™ includes a proprietary risk algorithm (patent pending) detecting conductive hearing loss with 94% accuracy. Users also have the option to add the [CEDRA](#) questionnaire for ear disease.



hearX Cloud


Effortlessly review patient results and provide remote support, counseling, assistance, and detailed reporting. Test data integrates seamlessly with hearX Cloud, our robust cloud-based data management system, ensuring streamlined workflows and enhanced efficiency.

Features you can benefit from




Environmental noise warning

Real-time monitoring of noise to alert the user of environmental noise concerns.



Adjustable protocols

Tailor your testing protocols to specific patient needs, ensuring a personalized and accurate hearing health assessment.



Centralized, real-time data management

Seamless integration into hearX Cloud, our cloud data management system.

And even more...

- Conditioning functionality
- Patient signature for consent
- Audiogram and degree of loss classification
- Narrow band masking





See our Self Test Kit in action

Schedule a risk-free demo now to learn how the hearX® Self Test Kit can make your hearing assessment processes more streamlined, cost-effective and easy to manage.


[Book a FREE demo now](#)

Technical specifications

- Pure tone test frequencies | 0.125 to 8 kHz
- Rise/fall time | ~35 ms
- Total harmonic distortion | <2%
- Test intensity levels | 10 to 90 dB HL

Take data management to the next level with hearX Cloud

hearX® Self Test Kit is managed through hearX Cloud, enabling seamless synchronization between hearTest on your hearX hardware set and hearX Cloud (Electronic Health Record Platform). This ensures that all test and patient records are accessible in a centralized online location.



www.hearxgroup.com | US: +1 (415) 212-5500 | ZA: +27 (0) 12 030-0268

OR email sales@hearxgroup.com

USA

2140 South DuPont Highway, Camden,
Kent County Camden, DE 19934


ZA

Ground Floor, Building 2, Ashlea Gardens
Office Park, 180 Garsfontein Road, Ashlea
Gardens, Pretoria, 0081

Please note: Shipping charges are included for local shipments in the USA. Additional charges will apply for shipments outside of the USA.

hearX® is committed to protecting your privacy and ensuring that personal information is collected and used lawfully and transparently. We use the information provided to contact you about relevant products, services and content. You may unsubscribe at any time. Standard hearX Privacy Policy applies.

FDA Registered | CE-marked | ISO 13485 Certified | SAHPRA Registered | HIPAA Compliant | POPIA Compliant | GDPR Compliant



Copyright © 2016 - 2025 hearX USA. Copyright in this work including its layout and the designs contained herein is reserved and the making of copies is prohibited.